


MR Elastography Phantom

Scanning Parameter Recommendations

- To set-up the phantom for scanning, refer to the illustrations in the MRE PHANTOM USER GUIDE included with each MRE Phantom.
- Use the following recommended MRE Phantom imaging parameters that are specific to your GE MR scanner.

| GE - 2D MRE Phantom Parameter Recommendations | | | |
|---|--|--|---|
| Scanners and Sequences (Note 1) | Scanner | Artist, Creator, Explorer, HDx, Optima MR450w, Voyager | Architect, Discovery MR750w, PET/MR, Pioneer, Premier |
| | Field Strength | 1.5T | 3T |
| | Software versions (Verified) | DV22.1, DV24, DV25, DV26 | DV22.1, DV24, DV25, |
| | Software versions (Compatible) | | |
| | Pulse sequence | MR-Touch (GRE) | MR-Touch (EPI) |
| | Mode | 2D | 2D |
| | Options | Fast, ASSET, MultiPhase | ASSET, FC |
| Phantom Setup | Phantom Setup | See MRE PHANTOM USER GUIDE with Phantom Kit. | |
| | Coil (Note 2) | Head | |
| Slice Positing | Place one coronal slice at the center of the height of the phantom | | |
| |  | | |
| Information Input (Pretend Patient) | Position | feet-first, supine | feet-first, supine |
| | Weight | 150 Lbs (68 kg) | 150 Lbs (68 kg) |
| | Height | 5 ft (1.5 m) | 5 ft (1.5 m) |
| Coil | Coil (Note 2) | head | head |
| Imaging Parameters | Imaging Plane | coronal | coronal |
| | Number of slices | 1 | 1 |
| | Slice thickness (mm)/gap | 10 mm / 0 mm | 8 mm / 2 mm |
| | FOV (mm) / Phase FOV (100%) (Note 3) | 20 cm/100% | 20 cm/100% |
| | Matrix | 256 × 64 | 32 × 32 |
| | TE (msec) | min full TE (typically ~18.2) | min full (around 57.6 msec) |

| | | | |
|------------------------------|--|---|---|
| Imaging Parameters | TR (msec) | 50 | 250 |
| | Flip Angle (degree) | 25 | 90 (default) |
| | NEX | 1 | 1 |
| | EPI shots | N/A | 1 |
| | Bandwidth (kHz) | 31.25 | 250 (hard coded) |
| | Freq Encoding Dir | Superior-Inferior | Superior-Inferior |
| | Phase Acq. Order | N/A | N/A |
| | Delay After Acq. | N/A | N/A |
| | Acceleration (Note 4) | ASSET | ASSET |
| | Acceleration factor (Note 4) | 1 | 2 |
| | Number of breath holds | N/A | N/A |
| | Shimming Volume | Cover the whole phantom | Cover the whole phantom |
| | Spectrum Peaks | Peak with middle freq (there are 3 peaks) | Peak with middle freq (there are 3 peaks) |
| | Saturation Band | SI | SI |
| | Scan Time | about 28 sec (Note 2) | 10 sec |
| MR-Touch Tab (Note 1) | Temporal Phases | 4 | 4 |
| | MEG Frequency (Hz) | 75 | 90 |
| | Driver Frequency (Hz) | 60 | 60 |
| | Driver Amplitude (%) (Note 6) | 10 | 10 |
| | Driver Cycle Per Trigger | 3 | 15 (Cannot edit) |
| | MEG Direction | Z | Z |
| | MENC um/rad | 28.5 typ (Cannot edit) | 28.5 typ (Cannot edit) |
| Advanced Tab | Use resoundant | 1 | 1 |
| | MRE Gradient Amplitude (g/cm) (Note 5) | 1.84 | 1.84 |

NOTES:

- (1) Specific tab and parameters vary based on different software versions and MRE sequences; the generic MRE parameters for driver and motion encoding gradients are the guideline to those specific tab and parameters (MRE-related); overall, this recommendation is conservative so that it can be successfully performed at all software versions and scanners.
- (2) Use of a multi-channel RX head coil is preferred. Alternatively, the Torso/Spine coils can be used.
- (3) 20cm FOV is a good minimum value for phantom studies, smaller FOVs may reduce SNR or confidence and should be avoided. Larger FOV may be used and may be beneficial for correlation to in vivo scans. FOV should be kept consistent across all phantom scans or results will not be comparable though time.
- (4) ASSET acceleration factor $R = 1$ is recommended for GRE, while $R = 2$ for EPI to allow shorter TE. $R > 1$ can be used generally, particularly for correlation to in vivo sequences. Larger FOVs may be desirable to avoid aliasing artifacts.
- (5) Performance may vary based on gradient hardware performance and maximum gradient amplitude. 20mT/m is a good value for 10% acoustic amplitude in phantom; 4mT/m is good for 50% acoustic amplitude.

Questions - Questions regarding the Resoundant MRE Phantom Scanning Parameter Recommendations may be directed to:

Resoundant, Inc.
 421 First Avenue SW STE 204W
 Rochester, Minnesota 55902 USA
 Phone: 507.322.0011
 Email: mreinfo@resoundant.com